

# PressureGuard® Bariatric APM

Category: Bariatric Alternating Pressure



## Feature-Packed Therapy for Users up to 750 lbs.

Segmented top reduces sharp pressure drop-offs, while air system below handles the “heavy lifting,” alternating weight-bearing across the surface in user-selectable intervals. “Two up, one down” inflation pattern ensures maximum patient support throughout the cycle.

### Powerful Digital Control Unit

- Up to four times the air output of conventional units.
- Eight electronically-controlled comfort settings.
- Programmable cycle times of 20, 25, 30 or 35 minutes.
- Low pressure indicator and audible alarm.



Additional bariatric options available – see page 22.

All models are fire barrier version (FBV) as indicated by “-29” suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.

### QUICK FACTS

## PressureGuard® Bariatric APM

### Primary Uses:

- Aggressive treatment (through Stage IV)
- Prevention of pressure ulcers in sedentary or high risk patients

### Patient Indications:

- Can’t or won’t reposition frequently
- Can be positioned off of pressure ulcer

### Key Advantages:

- Simple operation
- Gentle and effective pressure changes
- Customizable firmness and time settings
- PressureGuard design: stability, safety, anti-shearing

Item #	Description
AP8042-29	Bariatric APM, 80”L x 42”W x 7”H
AP8048-29	Bariatric APM, 80”L x 48”W x 7”H
AP8054-29	Bariatric APM, 80”L x 54”W x 7”H
* AP8442-29	Bariatric APM, 84”L x 42”W x 7”H
* AP8448-29	Bariatric APM, 84”L x 48”W x 7”H
* AP8454-29	Bariatric APM, 84”L x 54”W x 7”H
7500	Bariatric APM Control Unit Only

Medicare Approved: E0277 (Group II)  
 Weight limit: 750 lbs.  
 Warranty: 18 months (all components)  
 Mattress Weight: Approx. 25 lbs.  
 UL listed: Yes

\* Allow up to 10 working days lead time for this size.

