



# Span Surfaces for Recessed Deck Beds

Specialty models to accommodate the unique challenges of recessed deck beds.




Featuring Span's proprietary Safety Edge™ two-piece bolster system, adapted to overcome the design restrictions inherent in recessed decks.



Inner bolster "cradle" deflects when patient sits on or moves toward edge of bed . . .




Compressible section to accommodate multiple lengths and foot positions



. . . which helps direct patient's center of gravity safely back toward center of bed.

## Choose from among these leading Span surfaces


**Geo-Mattress® UltraMax™:** Premium therapeutic foam  
(Item # UMX85V-29  
– see page 24)



**PressureGuard® Custom Care®:** Non-powered reactive pressure redistribution support surface  
(Item # NP85V29  
– see page 16)



**PressureGuard® Span-Care® Convertible:** Non-powered reactive pressure redistribution support surface  
(Item # SC85V29  
– see page 13)



**PressureGuard® Custom Care® Convertible:** Convertible treatment with add-on alternating pressure and lateral rotation  
(Item # CJ85V29;  
add-on control unit  
item #6500  
– see page 14)



**PressureGuard® Custom Care® Convertible LAL:** Convertible treatment with add-on microclimate management, alternating pressure and lateral rotation  
(Item # CL85V29;  
add-on control unit  
item #8400  
– see page 15)



**PressureGuard® Easy Air® LR:** Microclimate management, alternating pressure and lateral rotation  
(Item # L85VLR-29  
– see page 5)

