

Fixed Ceiling Lifts

Fixed lifts for Healthcare Facilities

All the features of the portable lifts, plus

- Patented FastTrack systems for all applications.
- Facility-friendly design, including fray-resistant strap, waterproof handset, smart spreader bar, sensible emergency cord.
- Easy servicing without removal from the track.
- Enhanced safety features for caregivers and patients.



SPECIFICATIONS	
Product weight	Manual model: 17 lb Motorized model: 19 lb
Safe Working Load (SWL)	Fixed: 286, 440, or 600 lb
Noise level	Maximum 54 dBA
Lifting speed	2"/sec at 286 lb 1.8"/sec at 440 lb 1.6"/sec at 600 lb
Battery type	Lithium-ion 25.2V 2500 mAh
Lifts per charge	Approximately 80 cycles of 24" at 200 lb 70 cycles of 24" at 286 lb 45 cycles of 24" at 440 lb 30 cycles of 24" at 600 lb
Battery charging	Full capacity in approximately 2 hours
Life expectancy	10 years with only one battery replacement
Warranty	3-year warranty
Electrical Class	Medical electrical equipment Class 1
Protection Class	Type BF
Compliance	Meets international certification and compliance standards for medical devices

See our complete line of Silvalea slings at silvaleaslings.com

