

Span F450T

Standing Transfer Aid

- For users requiring Minimum Assistance
- 450 lb. (205 kg) safe working load
- Keeps users actively engaged in the transfer process
- **Easy base width adjustment** accommodates chairs and commodes.
- **A convenient toileting alternative** vs. use of a wheelchair.
- **Locking rear casters** add an extra measure of safety.
- **Attachment hooks** for optional support belt.
- **Simple, fast assembly** using tools provided.



SPECIFICATIONS

Weight Capacity	450 lb / 205 kg
Overall Height	43.2"
Overall Base Length	35.5"
Minimum Seat Height	26"
Maximum Seat Height	33.5"
Front-end Base Height	3.6"
Internal Leg Width (Open)	32"
External Leg Width (Closed)	27.7"
Horizontal Handlebar Width	20"
Net Weight	78 lb / 35.2 kg
Shipping Weight	93 lb / 42 kg

Item # Description

F450T	Sit-To-Stand Patient Lift
Warranty	3 years (frame)

Span F500P

Full Body Patient Lift

- For users requiring Maximum Assistance
- 500 lb. (228 kg) safe working load
- 6-point spreader bar for maximum sling options
- **Easy base width adjustment** with locked position; accommodates chairs and commodes.
- **Redundant operation** via control unit, allowing use if pendant is lost or damaged.
- **Floor extraction:** Can be used to lift patients from floor.
- **Overload protection** with visual/audible alert, halts lifting when capacity is exceeded.
- **Manual emergency lowering** for added safety.
 - **Quick-release battery** can be charged via charging cradle (optional), or via plug-in charger (included) for charging in place.
 - **Backlit LCD display** during operation for ease of viewing.
 - **Simple, fast assembly** using tools provided.



SPECIFICATIONS

Weight Capacity	500 lb / 228 kg
Base Legs Length	44.9"
Base Closed External Width	24.9"
Base Open External Width	42.9"
Mast Height	48.3"
Base Height	4.5"
Minimum Lifting Height	25.8"
Maximum Lifting Height	69.1"
Net Weight	99 lb / 45 kg
Shipping Weight	108 lb / 49 kg

Item # Description

F500P	Full Body Patient Lift
Warranty	3 years (frame and spreader bar); 2 years (electronics); 1 year (battery)

