



Span Surfaces for Recessed Deck Beds

Specialty models to accommodate the unique challenges of recessed deck beds.




Featuring Span's proprietary Safety Edge™ two-piece bolster system, adapted to overcome the design restrictions inherent in recessed decks.



Inner bolster "cradle" deflects when patient sits on or moves toward edge of bed . . .



Compressible section to accommodate multiple lengths and foot positions



. . . which helps direct patient's center of gravity safely back toward center of bed.

Choose from among these leading Span surfaces



PressureGuard® Easy Air® LR:
Microclimate management, alternating pressure and lateral rotation (Item # L85VLR-29 – see page 5)



Geo-Mattress® UltraMax™:
Premium therapeutic foam (Item # UMX85V-29 – see page 26)



PressureGuard® Custom Care® Convertible LAL:
Convertible treatment with add-on microclimate management, alternating pressure and lateral rotation (Item # CL85V29; add-on control unit item #8400 – see page 17)



PressureGuard® Custom Care® Convertible:
Convertible treatment with add-on alternating pressure and lateral rotation (Item # CJ85V29; add-on control unit item #6500 – see page 16)



PressureGuard® Custom Care®: Non-powered reactive pressure redistribution support surface (Item # NP85V29 – see page 18)

