



# Span Surfaces for Recessed Deck Beds

Specialty models to accommodate the unique challenges of recessed deck beds.




Featuring Span's proprietary Safety Edge™ two-piece bolster system, adapted to overcome the design restrictions inherent in recessed decks.



Inner bolster "cradle" deflects when patient sits on or moves toward edge of bed . . .



Compressible section to accommodate multiple lengths and foot positions



. . . which helps direct patient's center of gravity safely back toward center of bed.

## Choose from among these leading Span surfaces



**PressureGuard® Easy Air® LR:**  
Microclimate management, alternating pressure and lateral rotation (Item # L85VLR-29 – see page 5)



**Geo-Mattress® UltraMax™:**  
Premium therapeutic foam (Item # UMX85V-29 – see page 26)



**PressureGuard® Custom Care® Convertible LAL:**  
Convertible treatment with add-on microclimate management, alternating pressure and lateral rotation (Item # CL85V29; add-on control unit item #8400 – see page 17)



**PressureGuard® Custom Care® Convertible:**  
Convertible treatment with add-on alternating pressure and lateral rotation (Item # CJ85V29; add-on control unit item #6500 – see page 16)



**PressureGuard® Custom Care®:** Non-powered reactive pressure redistribution support surface (Item # NP85V29 – see page 18)

