


# Bariatric Solutions

## Four Options for Weights up to 1000 lbs.




**PressureGuard Easy Air® XL** low air loss/alternating pressure mattress system  
(see page 4 and 5 for series description)

Part#	Description	Weight range
L8042XL-29	Easy Air XL system, 80"L x 42"W x 7"H	250 - 1000 lbs.
L8048XL-29	Easy Air XL system, 80"L x 48"W x 7"H	250 - 1000 lbs.
* L8053XL-29	Easy Air XL system, 80"L x 53"W x 7"H	250 - 1000 lbs.
* L8442XL-29	Easy Air XL system, 84"L x 42"W x 7"H	250 - 1000 lbs.
* L8448XL-29	Easy Air XL system, 84"L x 48"W x 7"H	250 - 1000 lbs.
* L8453XL-29	Easy Air XL system, 84"L x 53"W x 7"H	250 - 1000 lbs.
8300	Easy Air XL Bariatric Control Unit Only	



**PressureGuard® Bariatric APM** alternating pressure air mattress system  
(see pages 10 and 11 for product description)

Part#	Description	Weight range
AP8042-29	Bariatric APM, 80"L x 42"W x 7"H	150 - 750 lbs.
AP8048-29	Bariatric APM, 80"L x 48"W x 7"H	150 - 750 lbs.
AP8054-29	Bariatric APM, 80"L x 54"W x 7"H	150 - 750 lbs.
* AP8442-29	Bariatric APM, 84"L x 42"W x 7"H	150 - 750 lbs.
* AP8448-29	Bariatric APM, 84"L x 48"W x 7"H	150 - 750 lbs.
* AP8454-29	Bariatric APM, 84"L x 54"W x 7"H	150 - 750 lbs.
7500	Bariatric APM Control Unit Only	



**PressureGuard CFT®** non-powered treatment mattress system  
(see page 12 for product description)

Part#	Description	Weight range
CF8042-29	PressureGuard CFT, 80"L x 42"W x 7"H	100 - 750 lbs.
CF8048-29	PressureGuard CFT, 80"L x 48"W x 7"H	100 - 750 lbs.
CF8054-29	PressureGuard CFT, 80"L x 54"W x 7"H	100 - 750 lbs.
* CF8442-29	PressureGuard CFT, 84"L x 42"W x 7"H	100 - 750 lbs.
* CF8448-29	PressureGuard CFT, 84"L x 48"W x 7"H	100 - 750 lbs.
* CF8454-29	PressureGuard CFT, 84"L x 54"W x 7"H	100 - 750 lbs.



**Geo-Mattress Atlas®** therapeutic foam mattress  
(see page 29 for product description)

Part#	Description	Weight limit
A8035-29	Geo-Mattress Atlas, 80"L x 35"W x 7"H	600 lbs.
A8042-29	Geo-Mattress Atlas, 80"L x 42"W x 7"H	750 lbs.
A8048-29	Geo-Mattress Atlas, 80"L x 48"W x 7"H	750 lbs.
A8054-29	Geo-Mattress Atlas, 80"L x 54"W x 7"H	750 lbs.
* A8442-29	Geo-Mattress Atlas, 84"L x 42"W x 7"H	750 lbs.
* A8448-29	Geo-Mattress Atlas, 84"L x 48"W x 7"H	750 lbs.
* A8454-29	Geo-Mattress Atlas, 84"L x 54"W x 7"H	750 lbs.

**Non-prorated warranty:** 2 years (Easy Air), 18 months (APM); 5 years (CFT), 3 years (Atlas).  
**Cover cleaning:** Wipe clean with gentle suds or approved cleaning agents. Avoid bleach.  
**\* Lead times:** Allow up to 10 working days from time of order for shipment for items noted by an asterisk (\*).  
**Shipping:** All models ship fully assembled in a flat carton via common carrier.  
**UL Listing:** APM conforms to UL 60601-1, CAN/CSA C22.2 NO.601.1.  
 Easy Air conforms to UL 60601-1, CAN/CSA C22.2 NO.601.1.

