



better science. better care.

PressureGuard® Air Therapy Surfaces Matrix

		Non-powered Air Therapy				Alternating/Lateral Rotation					Low Air Loss			Bariatric			
		Custom Care®	CFT®	Custom Care® Convertible	Custom Care® Convertible LAL	APM	APM2™	Safety Supreme™	APM3™	Custom Care® Convertible with control unit	Turn Select®	Easy Air®	Easy Air® LR	Custom Care® Convertible LAL with control unit	Bariatric CFT®	Bariatric APM	Easy Air® XL
Therapy Modalities	Non-powered pressure redistribution	✓	✓	✓	✓										✓		
	Convertible to powered via add-on control unit			✓	✓												
	Alternating Pressure					✓	✓	✓	✓	✓		✓	✓	✓		✓	✓
	Lateral Rotation (P) = Programmable						✓	✓	✓	✓	✓ (P)		✓	✓			
	Low Air Loss											✓	✓	✓			✓
	Powered Flotation										✓	✓	✓	✓			✓
Frame Types	Versions to accommodate CareAssist® and recessed deck frames	✓		✓	✓					✓			✓				
Control Unit Features	Alarm: (V)isual, (A)udible			V, A	V		V, A	V, A	V, A	V, A	V, A	V	V	V		V, A	V
	Timed Auto Firm (20 minutes)							✓	✓								
Microclimate Management	(P) = Passive, (A) = Active	P	P	P	P	P	P	P	P	P	A	A	A	P	P	A	
Mattress Construction	Geo-Matt® top surface	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Safety Edge™ bolster design (RP) = Raised perimeter	✓	✓	✓	✓	✓	✓	✓ (RP)	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Can maintain air in power outage; (A) = Automatically	✓(A)	✓(A)	✓(A)	✓(A)	✓	✓	✓	✓(A)	✓(A)	✓(A)	✓(A)	✓(A)	✓(A)	✓	✓	✓(A)
	Heel Slope™; Foam-only heel section	✓(F)	✓	✓(F)	✓(F)				✓(F)	✓(F)				✓(F)			
	100% Latex free, 100% DHEP-free	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Cover	Non-stretch nylon with butyl/neoprene coating		✓			✓	✓	✓			✓				✓	✓	
	Bi-directional stretch with Shear Transfer Zone™ design	✓		✓					✓	✓							
	Zippered air flow coverlet with Air Diffusion Matrix™ layer				✓							✓	✓	✓			✓
General	Medicare Code		E0373	E0373	E0373	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0373	E0277	E0277
	Therapeutic Weight Limit (lbs.)	500	500	500	500	350	350	350	350	350	350	500	500 (350 LR)	500 (350 LR)	750	750	1000
	Warranty: Mattress, non-pro-rated	5 yrs.	5 yrs.	5 yrs.	5 yrs.	18 mos.	18 mos.	18 mos.	2 yrs.	5 yrs.	2 yrs.	2 yrs.	2 yrs.	5 yrs.	5 yrs.	18 mos.	2 yrs.
	Warranty: Control unit, non-pro-rated					18 mos.	18 mos.	18 mos.	2 yrs.	2 yrs.	2 yrs.	2 yrs.	2 yrs.	2 yrs.		18 mos.	2 yrs.
	Meets Fire Safety standards Cal 129, CPSC16CFR1633, Boston IX-11	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

*Hill-Rom is a mark, and CareAssist® is a registered mark, of Hill-Rom® Services, Inc.