

Seating

Category: Prevention without Positioning

Gel-T®

Segmented gel-filled bladder cradled between two soft foam layers. Provides excellent pressure redistribution and pain relief for those at moderate risk of pressure ulcers. Wipe-clean cover keeps out moisture and helps reduce shear and friction.

Cover also features non-skid bottom. Perfect for use in home or facility; 18-month warranty. Height 2.5".

Weight limit 250 lbs.

* sizes 18" x 18" and larger up to 350 lbs. (500 lbs. on Bariatric sizes)

Medicare Code:

General Use – E2601.

For sizes 22" or wider, use E2602.



| Item # | Description | Qty/Cs |
|------------------|---------------------------------------|--------|
| 8015-05 | Gel-T, 15"W x 15"L | 1 |
| 8016-05 | Gel-T, 16"W x 16"L | 1 |
| 801618-05 | Gel-T, 16"W x 18"L | 1 |
| 8017-05 | Gel-T, 17"W x 17"L | 1 |
| 801816-05 | Gel-T, 18"W x 16"L | 1 |
| 8018-05 | Gel-T, 18"W x 18"L | 1 |
| 8019-05 | Gel-T, 19"W x 19"L (fits Geri-chairs) | 1 |
| 802018-05 | Gel-T, 20"W x 18"L (fits Geri-chairs) | 1 |
| 8021-05 | Gel-T, 21"W x 21"L (fits Geri-chairs) | 1 |
| Bariatric sizes: | | |
| B802218-05 | Gel-T, 22"W x 18"L | 1 |
| B802418-05 | Gel-T, 24"W x 18"L | 1 |

Fire barrier versions can be ordered to conform with 2000 NFPA 101 (Life Safety Code) and Cal TB# 133. To order, replace the "-05" suffix with "-29" in item

QUICK FACTS

Use for:

- Comfort
- Prevention of breakdown

RISK LEVEL



Geo-Matt® PRT (3" high density)

Modification of Geo-Matt's time-tested design creates lower, firmer base needed by more active wheelchair users. Same individually responsive cells, Ring-of-Air® heat/moisture reduction, and pressure redistribution as original Geo-Matt. Cover consists of anti-slip top surface and non-skid bottom. Built-in sleeve in bottom of cover allows use of molded Sling-Fill base (not included-see page 33); 18-month warranty. Height 3".

Weight limit 250 lbs.

(4" height/500 lb. weight limit on Bariatric sizes)

Medicare Code:

Skin Protection – E2603

(only versions with fabric cover).

For sizes 22" or wider, use E2604.



| Item # | Description | Qty/Cs |
|----------------------------|---|--------|
| 50746-070 | Geo-Matt PRT, 18"W x 16"L, w/ Soft Skin® sleeve | 6 |
| 50746-083 | Geo-Matt PRT, 18"W x 16"L, w/ vinyl/nylon cover | 6 |
| PRT1616-96 | Geo-Matt PRT, 16"W x 16"L, w/ anti-slip cover | 6 |
| PRT1616-96-1 | Geo-Matt PRT, 16"W x 16"L, w/ anti-slip cover | 1 |
| PRT1816-96 | Geo-Matt PRT, 18"W x 16"L, w/ anti-slip cover | 6 |
| PRT1816-96-1 | Geo-Matt PRT, 18"W x 16"L, w/ anti-slip cover | 1 |
| PRT1818-96 | Geo-Matt PRT, 18"W x 18"L, w/ anti-slip cover | 6 |
| PRT1818-96-1 | Geo-Matt PRT, 18"W x 18"L, w/ anti-slip cover | 1 |
| PRT2018-96 | Geo-Matt PRT, 20"W x 18"L, w/ anti-slip cover | 6 |
| PRT2018-96-1 | Geo-Matt PRT, 20"W x 18"L, w/ anti-slip cover | 1 |
| Bariatric sizes (4" high): | | |
| BPRT2218-96 | Geo-Matt PRT, 22"W x 18"L, w/anti-slip cover | 6 |
| BPRT2218-96-1 | Geo-Matt PRT, 22"W x 18"L w/anti-slip cover | 1 |
| BPRT2418-96 | Geo-Matt PRT, 24"W x 18"L w/anti-slip cover | 6 |
| BPRT2418-96-1 | Geo-Matt PRT, 24"W x 18"L w/anti-slip cover | 1 |

Fire barrier versions can be ordered to conform with 2000 NFPA 101 (Life Safety Code) and Cal TB# 133. To order, replace the "-96" suffix with "-29" in item number.

QUICK FACTS

Use for:

- Comfort
- Prevention of breakdown

RISK LEVEL



Geo-Matt® (4" medium density)

Proven design provides outstanding pressure redistribution and comfort. Unique geometry creates individually responsive cells that reduce shearing and cradle bony prominences. Ring-of-Air® design helps reduce heat and moisture build-up; 18-month warranty. Height 4".

All models 18"W x 16"L.

Weight limit 200 lbs.

Medicare Code:

Skin Protection – E2603 (only versions with fabric cover).



| Item # | Description | Qty/Cs |
|-----------|--|--------|
| 50746-040 | Geo-Matt, cushion only | 6 |
| 50746-050 | Geo-Matt with Soft Skin® sleeve | 6 |
| 50746-051 | Geo-Matt with fire retardant Soft Skin® sleeve | 6 |
| 50746-084 | Geo-Matt with anti-slip cover | 6 |

QUICK FACTS

Use for:

- Comfort
- Prevention of breakdown

RISK LEVEL

