


Encore Structural Highlights



Encore™

Optional Accessories:

- Bed extenders to 84"
- Ergonomic 3-position assist device
- Bed & deck angle indicator
- Snap-in pendant holder fits all rails
- Molded plastic half-rail
- Headboard mount pendant holder
- Transport dollies
- Portable backup battery
- Regular backup battery
- Drainage hooks
- IV pole
- Trapeze

Height Floor to Mattress Support	7-7/8" - 31"
Deck Angles - Back Elevation	0 - 70°
Deck Angles - Knee Elevation	0 - 30°
Bed Length - Overall Length	82", 86", 90"
Bed Length - Retraction Travel	12"
Bed Length - Sleep Surface	76", 80", 84"
Bed Width - With Side Rails	40-1/2", 43-1/2", 46-1/2"
Bed Width - Sleep Surface	35", 39", 42"
Weight of Bed (without accessories)	253 lbs / 115 kg
Safe Working Load	600 lbs / 272 kg

Item #	Description
MC9R	Encore bed with GlideAlign retractability
MC9RW	Encore bed with GlideAlign retractability, with readywide expanders up to 42" wide
MC9RZ	Encore bed with GlideAlign retractability, with built-in battery backup
MC9RWZ	Encore bed with GlideAlign retractability, with readywide expanders up to 42" wide and built-in battery backup
R9938	Three position pivot assist - universal
R9939	Three position pivot assist - set
R9308H	Half head end side assist/rail - set
R9308H-L	Half head end side assist/rail - left
R9308H-R	Half head end side assist/rail - right
Q6690	Extender kit 4" - to 84"
Q6699	Dolly set
Q6730	Universal Trapeze bar
Q6731	Encore Trapeze mounting frame
Q6708	Roller bumper kit
Q6703	Deck angle indicators with brackets



For ordering and technical support on beds, accessories and furniture, call: **800-268-8671**