



better science. better care.

## Seating Matrix Innovative Solutions

		Treatment w/Positioning		Positioning	Skin Protection			General Use							
															
		Isch-Dish®	Equalizer®	Isch-Dish® Thin	EZ Dish®	Geo-Matt® Contour	Geo-Matt® PRT	Geo-Matt®	Geo-Matt® Wedge	Gel-T®	Soft Pedal	Cross Cut	Short Wave®	Geo-Wave®	
Level of Care	General Use					✓			✓	✓	✓	✓	✓	✓	
	Skin Protection	✓	✓		✓	✓	✓	✓							
	Positioning Features	✓	✓	✓	✓	✓			✓						
	Treatment Surface	✓	✓												
Core	Geo-Matt®					✓	✓	✓	✓				✓	✓	
	Passive Microclimate Mgmt.					✓	✓	✓	✓				✓	✓	
	Minimizes Shearing					✓	✓	✓	✓			✓	✓	✓	
	Gel Components		✓							✓	✓				
Cover	Anti-Slip Bottom				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Machine Washable	✓		✓											
	Two-way Stretch Fabric Top	✓	✓	✓											
	Wipe Clean				✓	✓	✓	✓	✓	✓	✓		✓	✓	
General	Medicare Code *Bariatric Code	E2607	E2607	E2603	E2603	E2603	E2603 E2604*	E2603	E2601	E2601 E2602*		E2601	E2601 E2611*		
	Therapeutic Weight Limit (lbs.)	250	250	250	250	250	250	200	250	250	250	250	300	300	
	Bariatric Weight Limit (lbs.)		500					500		500					
	Warranty	2 yrs.	2 yrs.	2 yrs.	18 mos.	18 mos.	18 mos.	18 mos.	18 mos.	18 mos.	18 mos.	18 mos.	1 yr.	18 mos.	18 mos.
	Seat and Back Cushion												✓	✓	
	Height of Cushion (Inches)	3.5"	3.5"	2.25"	3"	3"	3" standard 4" bariatric	4"	4" front 2" back	2.5"	1.9"	3"	3"	3"	3"
	Viscoelastic Top		✓		✓						✓				