


Tool-less 42" Expander Kits for Rexx and Advantage Beds

Easy on, easy off expander kits install instantly. Require no tools, bolts, or pins. Each kit contains 6 deck expanders (3 each right and left) required to expand a 35" bed to 42" wide.

Bed frame, overhead view



Convenient decals indicate proper placement of each extender.



Simply snap the deck expander into place on both sides of the head, seat and knee decks (supplied in two sets of three).



A powerful latch beneath each expander unit firmly attaches to the deck making deck expansion to 42" a snap!



Molded, half-head assist rails can be mounted securely to the expanders before installation.



Kit is supplied with protective packaging for storage and re-use.



Item #	Description
Q6709	42" Deck Expander Kit for Rexx, Fast Rexx, and Advantage models

Item #	Description
Q6729	Universal Equipment Bar (foot end)



For ordering and technical support on beds, accessories and furniture, call: **800-268-8671**