

Span+Aids® Patient Positioners

Category: Extremities



Foot Drop Stop – Standard

Elevates heels, supports ankle position. Adjustable strap can accommodate ankle contracture, protect against foot drop. Square base provides stability. Ring-of-Air® design provides ventilation.
50638-511 *6/cs



Limb Elevator

Elevates casted or uncasted leg. Can be used in both supine and side lying positions.
50632-508 *6/cs
50632-308 2/cs




Foot Drop Stop – Original

Similar to standard version (above), except hexagonal base maximizes comfort by allowing user to lie with leg in multiple positions.
55555-01 *6/cs
50638-311 6/cs




Flared Leg Extension

Elevates leg and knee in extended position.
SP291-000 2/cs



Cradle Boot – Original

Floats heel, supporting foot in neutral position, prevents external leg rotation.
55555-02 *6/cs



Leg Trough

Immobilizes leg and knee in neutral position. Base is contoured to posterior leg.
SP033-000 2/cs



Cradle Boot – Standard

Floats heel, supporting foot in the neutral position. Well-ventilated boot is easy to use and apply. Square base provides stability, discouraging leg rotation.
50644-515 *6/cs



Upper Limb Positioner

Positions the entire upper arm with sloped groove for stability and comfort. Initially designed for use with burn patients.
SP121-000 2/cs



Cock-Up Splint

Designed to help maintain hand and wrist in a desired position. Holds fingers and thumb opposed to prevent contractures. Hand stirrup is adjustable.
50640-512 *8/cs