

**SPAN**  
better science. better care.  
**PressureGuard®  
Air Therapy  
Surfaces Matrix**

	Non-powered Air Therapy				Alternating/Lateral Rotation				Microclimate Management ["Low Air Loss" (LAL)]				Bariatric			
	Custom Care®	CFI®	SpanCare® Convertible	Custom Care® Convertible	APM	APM™ / Safety Supreme	SpanCare® Convertible w/control unit	Custom Care® Convertible w/control unit	Protocol®	Easy Air® Easy Air® LR	Custom Care® Convertible LAL w/control unit	Bariatric CFI®	Bariatric APM	Easy Air® XL		
Therapy Modalities	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Non-powered pressure redistribution	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Convertible to powered via add-on control unit	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Alternating Pressure					✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Lateral Rotation (P) = Programmable					✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Microclimate Management (LAL)					✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Powered Flotation					✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Frame Types	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Versions to accommodate recessed deck frames	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Control Unit Features				V, A	V	V, A	V, A	V, A	V, A	V	V	V, A	V, A	V		
Alarm: (V)visual, (A)audible				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Timed Auto Firm (20 minutes)				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Microclimate Management	P	P	P	P	P	P	P	A	A	A	A	P	P	A		
(P) = Passive, (A) = Active																
Geo-Matt® top surface	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Safety Edge™ bolster design (RP) = Raised perimeter	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Mattress Construction	✓(A)	✓(A)	✓(A)	✓(A)	✓	✓	✓(A)	✓	✓	✓(A)	✓(A)	✓	✓	✓(A)		
Can maintain air in power outage; (A) = Automatically																
Heel Slope™; (F) = Foam-only heel section	✓(F)	✓	✓(F)	✓(F)	✓	✓	✓(F)	✓	✓	✓	✓(F)	✓	✓	✓		
100% Latex free, 100% DHEP-free	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Non-stretch nylon with butyl/neoprene coating	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Bi-directional stretch with Shear Transfer Zone™ design	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Cover																
Zippered air flow coverlet with Air Diffusion Matrix™ layer				✓				✓	✓	✓	✓	✓	✓	✓		
Medicare Code	E0373	E0373	E0373	E0373	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0373	E0277	E0277		
Therapeutic Weight Limit (lbs.)	500	500	500	500	350	350	350	400	500 (350 LR)	500 (350 LR)	500 (350 LR)	750	750	1000		
Warranty: Mattress, non-pro-rated	5 yrs.	5 yrs.	5 yrs.	5 yrs.	18 mos.	18 mos.	18 mos.	18 mos.	18 mos.	2 yrs.	5 yrs.	18 mos.	18 mos.	2 yrs.		
Warranty: Control unit, non-pro-rated					18 mos.	18 mos.	18 mos.	18 mos.	18 mos.	2 yrs.	2 yrs.	18 mos.	18 mos.	2 yrs.		
Meets Federal Fire Safety Regulations	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		

\*Meets coding requirements. Code applied for.