

	Non-powered Air Therapy					Alternating/Lateral Rotation				Microclimate Management ["Low Air Loss" (LAL)]				Bariatric		
	Custom Care®	CFT®	SpanCare® Convertible	Custom Care® Convertible	Custom Care® Convertible LAL	APM	APM™ / Safety Supreme	SpanCare® Convertible w/control unit	Custom Care® Convertible w/control unit	Protocol®	Easy Air® EasyAir® LR	Custom Care® Convertible LAL w/control unit	Bariatric CFT®	Bariatric APM	Easy Air® XL	
Therapy Modalities	✓	✓	✓	✓	✓								✓			
Non-powered pressure redistribution	✓	✓	✓	✓	✓											
Convertible to powered via add-on control unit			✓	✓	✓											
Alternating Pressure						✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Lateral Rotation (P) = Programmable						✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Microclimate Management (LAL)																
Powered Flotation						✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Versions to accommodate recessed deck frames	✓		✓	✓	✓											
Alarm: (V)isual, (A)udible				V, A	V	V, A	V, A	V, A	V, A	V, A	V	V	V, A	V, A	V	
Timed Auto Firm (20 minutes)						✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Microclimate Management (P) = Passive, (A) = Active	P	P	P	P	P	P	P	P	P	A	A	A	P	P	A	
Geo-Matt® top surface	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Safety Edge™ bolster design (RP) = Raised perimeter	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Can maintain air in power outage; (A) = Automatically	✓(A)	✓(A)	✓(A)	✓(A)	✓(A)	✓	✓	✓	✓(A)	✓	✓(A)	✓	✓	✓	✓(A)	
Heel Slope™; (F) = Foam-only heel section	✓(F)	✓	✓(F)	✓(F)	✓(F)	✓	✓	✓	✓(F)	✓	✓(F)	✓	✓	✓	✓	
100% Latex free, 100% DHEP-free	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Non-stretch nylon with butyl/neoprene coating	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Bi-directional stretch with Shear Transfer Zone™ design	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Zippered air flow coverlet with Air Diffusion Matrix™ layer																
Medicare Code	E0373	E0373	E0373*	E0373	E0373	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0373	E0277	E0277	
Therapeutic Weight Limit (lbs.)	500	500	500	500	500	350	350	350	350	400	500 (350 LR)	500 (350 LR)	750	750	1000	
Warranty: Mattress, non-pro-rated	5 yrs.	5 yrs.	5 yrs.	5 yrs.	5 yrs.	18 mos.	18 mos.	18 mos.	18 mos.	18 mos.	2 yrs.	2 yrs.	5 yrs.	18 mos.	2 yrs.	
Warranty: Control unit, non-pro-rated						18 mos.	18 mos.	18 mos.	18 mos.	18 mos.	2 yrs.	2 yrs.	18 mos.	18 mos.	2 yrs.	
Meets Federal Fire Safety Regulations	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	

*Meets coding requirements. Code applied for.