


# PressureGuard Easy Air®

Category: Low Air Loss with Alternating Pressure

## Superior Microclimate Management, None of the Compromises.



Easy Air's patented "Air Diffusion Matrix" design maximizes removal of excess moisture.



Clinicians choose low air loss systems to address maceration, the weakening of the skin due to excess moisture. The PressureGuard Easy Air is independently documented to outperform every leading low air loss system in this key performance area. Remarkably, it does so with none of the compromises low air loss prescribers have come to accept, thanks to its exclusive Air Diffusion Matrix™ technology and patented PressureGuard air/foam design.

Easy Air maintains inflation in a power outage. Its stable surface won't "swallow up" the user. And its Safety Edge™ won't collapse, making transfers safer. Doesn't make the user too hot or too cold, either. Ultra-simple to set up, clean and use. Selectable comfort level. "On/Off" buttons control air flow and alternating pressure features separately. System includes mattress, controller unit and machine washable coverlet.

### QUICK FACTS

## PressureGuard Easy Air®

#### Primary Use:

- Aggressive Treatment (through Stage IV)

#### Patient Indications:

- Excessive maceration contributes to development of skin breakdown
- Any stage ulcer where maceration is a contributing factor
- Prevention for those at high risk of skin breakdown (e.g. Braden score 12 or less)

#### Key Advantages:

- Documented superiority in moisture removal
- Won't deflate in power outage
- Includes alternating pressure mode
- PressureGuard® design: stability, safety, anti-shearing

Item #	Description
L7535-29	Easy Air System 75"L x 35"W x 7"H
* L7835-29	Easy Air System 78"L x 35"W x 7"H
L8035-29	Easy Air System 80"L x 35"W x 7"H
L8435-29	Easy Air System 84"L x 35"W x 7"H
* L8032-29	Easy Air System 80"L x 32"W x 7"H
* L8039-29	Easy Air System 80"L x 39"W x 7"H
8000	Easy Air Controller Unit
CLT-L7535	Easy Air Coverlet 75"L x 35"W
* CLT-L7835	Easy Air Coverlet 78"L x 35"W
CLT-L8035	Easy Air Coverlet 80"L x 35"W
CLT-L8435	Easy Air Coverlet 84"L x 35"W
* CLT-L8032	Easy Air Coverlet 80"L x 32"W
* CLT-L8039	Easy Air Coverlet 80"L x 39"W


\* Allow up to 10 working days lead time for this size.

#### Complete System Includes Mattress, Control Unit, and Coverlet.

To order Coverlet only, add "CLT-" to part number (example: CLT-L8035).

Medicare Approved:	E0277 (Group II)
Weight limit:	500 lbs.
Warranty:	2 years, non-prorated
Mattress Weight:	22 lbs.
UL listed:	Yes

All models are fire barrier version (FBV) as indicated by "-29" suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.



Data Source:  
A Proposed Method for Quantifying Low-Air-Loss Mattress Performance by Moisture Transport.  
Richard S. Figliola PhD, PE, Clemson University, Clemson, SC, *Ostomy/Wound Management*, 2003;49(1):32-41.

