

Therapy Provided

- ✓ Microclimate Management (LAL)
- ✓ Alternating Pressure
- ✓ Lateral Rotation
- ✓ Powered Pressure Redistribution


QUICK FACTS

Indications

(in addition to those for Easy Air, p. 4):

Lateral Rotation for patients who find frequent manual repositioning painful or need more frequent repositioning

PressureGuard EASY AIR® LR



The all-in-one Easy Air LR model puts caregiver just a button-push away from a powerful, additional treatment modality, "Lateral Rotation" (LR). The unique PressureGuard approach to LR has been documented to help wounds close more than three times faster than on other commonly provided air therapy surfaces¹. Beyond pure pressure management, Lateral Rotation gives caregivers a versatile tool that is indicated for a wide variety of user needs.

1- Anderson, C. et al. Lateral rotation mattresses for wound healing. *Ostomy/Wound Management* April 2004;50(4):50-62.

Item #	Description
L7535LR-29	Easy Air LR System 75"L x 35"W x 7"H
L8035LR-29	Easy Air LR System 80"L x 35"W x 7"H
L8435LR-29	Easy Air LR System 84"L x 35"W x 7"H
L85VLR-29	Easy Air LR for recessed deck frame (see page 19)
* L8039LR-29	Easy Air LR System 80"L x 39"W x 7"H
* L8053LR-29 8200	Easy Air LR System 80"L x 53"W x 7"H Easy Air LR control unit only

*Allow up to 10 working days lead time for this size.

Medicare Approved: E0277 (Group II)
Weight limit: LR models support 500 lbs in FLOAT and ALTERNATING, 350 lbs in ROTATION.
Warranty: 2 years, (all components)
Mattress Weight: Approx. 22 lbs.
UL listed: Yes

PressureGuard EASY AIR® XL



Therapy Provided

- ✓ Microclimate Management (LAL)
- ✓ Alternating Pressure
- ✓ Powered Pressure Redistribution

QUICK FACTS

Indications:

(in addition to those for Easy Air, p. 4):

Patients weighing 500-1000 pounds with excess skin moisture

Item #	Description
L8042XL-29	Easy Air XL system, 80"L x 42"W x 7"H
L8048XL-29	Easy Air XL system, 80"L x 48"W x 7"H
* L8053XL-29	Easy Air XL System, 80"L x 53"W x 7"H
* L8442XL-29	Easy Air XL System, 84"L x 42"W x 7"H
* L8448XL-29	Easy Air XL System, 84"L x 48"W x 7"H
* L8453XL-29 8300	Easy Air XL System, 84"L x 53"W x 7"H Easy Air XL Bariatric control unit only

* Allow up to 10 working days lead time for this size.

All models are fire barrier version (FBV) as indicated by "-29" suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.

Additional bariatric options available – see page 20.

