


PressureGuard Easy Air® LR

Category: Low Air Loss with Alternating Pressure and Lateral Rotation

Unmatched Microclimate Management – Plus Two Dynamic Therapies.

The all-in-one Easy Air LR model puts caregiver just a button-push away from a powerful, additional treatment modality: Continuous Lateral Rotation Therapy (CLRT). The unique PressureGuard approach to CLRT has been documented to help wounds close more than three times faster than on other commonly provided air therapy surfaces¹. But beyond pure pressure management, CLRT gives caregivers a versatile tool that is indicated for a wide variety of user needs.

1- Anderson, C. et al. Lateral rotation mattresses for wound healing. Ostomy/Wound Management April 2004;50(4):50-62.



Item #	Description
L7535LR-29	Easy Air LR System 75"L x 35"W x 7"H
L8035LR-29	Easy Air LR System 80"L x 35"W x 7"H
L8435LR-29	Easy Air LR System 84"L x 35"W x 7"H
* L8039LR-29	Easy Air LR System 80"L x 39"W x 7"H
* L8053LR-29	Easy Air LR System 80"L x 53"W x 7"H
8200	Easy Air LR Control Unit
Medicare Approved:	E0277 (Group II)
Weight limit:	LR models support 500 lbs in FLOAT and ALTERNATING, 350 lbs in ROTATION.
Warranty:	2 years, (all components)
Mattress Weight:	Approx. 22 lbs.
UL listed:	Yes

PressureGuard Easy Air® XL

Category: Bariatric Low Air Loss up to 1000 lbs.



Heavy Duty micro-climate management with alternating pressure. Minimizes excess skin moisture – especially important for at-risk obese patients.

Item #	Description
L8042XL-29	Easy Air XL system, 80"L x 42"W x 7"H
L8048XL-29	Easy Air XL system, 80"L x 48"W x 7"H
* L8053XL-29	Easy Air XL System, 80"L x 53"W x 7"H
* L8442XL-29	Easy Air XL System, 84"L x 42"W x 7"H
* L8448XL-29	Easy Air XL System, 84"L x 48"W x 7"H
* L8453XL-29	Easy Air XL System, 84"L x 53"W x 7"H
8300	Easy Air XL Bariatric Control Unit Only

All models are fire barrier version (FBV) as indicated by "-29" suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.

Additional bariatric options available – see page 22.

