

Bariatric Solutions

Four Options for Weights up to 1000 lbs.



PressureGuard Easy Air® XL low air loss/alternating pressure mattress system (see page 5 for product description)

Part#	Description	Weight limit
L8042XL-29	Easy Air XL system, 80"L x 42"W x 7"H	1000 lbs.
L8048XL-29	Easy Air XL system, 80"L x 48"W x 7"H	1000 lbs.
8300	Easy Air XL Bariatric Control Unit Only	

Optional Pendant controller:
P02809 Wired, hand-held pendant for XL models only




PressureGuard® APM Bariatric alternating pressure air mattress system (see page 7 for product description)

Part#	Description	Weight limit
AP8042-29	PressureGuard APM, 80"L x 42"W	750 lbs.
AP8048-29	PressureGuard APM, 80"L x 48"W	750 lbs.
AP8054-29	PressureGuard APM, 80"L x 54"W	750 lbs.
7500	Bariatric APM Control Unit Only	

Powerful Digital Control Unit



- Up to four times the air output of conventional units
- Eight electronically-controlled comfort settings
- Programmable cycle times of 20, 25, 30 or 35 minutes
- Low pressure indicator and audible alarm



PressureGuard CFT® non-powered treatment mattress system (see page 8 for product description)

Part#	Description	Weight limit
CF8042-29	PressureGuard CFT, 80"L x 42"W	750 lbs.
CF8048-29	PressureGuard CFT, 80"L x 48"W	750 lbs.
CF8054-29	PressureGuard CFT, 80"L x 54"W	750 lbs.



Geo-Mattress Atlas® therapeutic foam mattress (see page 12 for product description)

Part#	Description	Weight limit
A8035-29	Geo-Mattress Atlas, 80"L x 35"W	600 lbs.
A8042-29	Geo-Mattress Atlas, 80"L x 42"W	750 lbs.
A8048-29	Geo-Mattress Atlas, 80"L x 48"W	750 lbs.
A8054-29	Geo-Mattress Atlas, 80"L x 54"W	750 lbs.

Non-prorated warranty: 18 months (APM); 3 years (Atlas), 5 years (CFT).
Cover cleaning: Wipe clean with gentle suds or approved cleaning agents. Avoid bleach.
Lead times: Allow up to 10 working days from time of order for shipment.
Shipping: All models ship fully assembled in a flat carton via common carrier.
UL Listing: APM conforms to UL 60602-1, CAN/CSA C22.2 NO.601.1.
 Easy Air conforms to UL 60601-1, CAN/CSA C22.2 NO.601.1.