



Geo-Mattress® UltraMax™



Now made with



bleach-resistant fabric cover
(see page 15 for details)



Designed to exceed even the most rigorous requirements of the acute care environment.

Precision-engineered to deliver critical performance advantages inside and out! On the outside, the individually articulating cells of the clinically proven Geo-Matt® top surface combine with Span-America's proprietary Shear Transfer Zones® design to provide unmatched protection from the harmful effects of micro shear, macro shear and rotational shear. Helps reduce caregiver injury by "gliding" the user back to position following head of bed elevation. On the inside, the 3-Dimensional Zoning™ core design of the UltraMax provides anatomically calibrated levels of support and unmatched user comfort. And in the seat area – where the mechanical fatigue of bed gatching and concentrated user weight affect durability most – the unique cell structure of the high resiliency UHP 500™ layer provides a 55% higher Support Factor than found in typical medical grade foams. Viscoelastic foam and gentle Heel Slope® provide added protection at vulnerable heels.



SURFACES

Therapy Provided

- ✓ Non-powered Pressure Redistribution

QUICK FACTS

Indications:

- Prevention and early intervention of pressure ulcers
- Moderate to high risk for skin breakdown
- Decreased mobility
- Repositions frequently by self or caregiver
- Treatment of stage 1-4 pressure ulcers¹

Key Advantages:

- Minimized shearing thru Geo-Matt® surface geometry and Shear Transfer Zones® within stretch cover (see description on page 15)
- Outstanding comfort for longer, more restful sleep
- Superior in-bed mobility, caregiver ergonomics
- Enhanced durability, protection from "bottoming out"

Item #	Description
UMX7536-29	75"L x 36"W x 6"H
UMX8036-29	80"L x 36"W x 6"H
UMX8436-29	84"L x 36"W x 6"H
UMX85V-29	For recessed deck frames (see page 17)
UMX80CA-29	80" L for CareAssist®+
UMX84CA-29	84" L for CareAssist®+
* UMX8039-29	80"L x 39"W x 6"H
* UMX8042-29	80"L x 42"W x 6"H
* UMX8439-29	84"L x 39"W x 6"H
* UMX8442-29	84"L x 42"W x 6"H

* Allow up to 10 working days lead time for these sizes.

Weight: 24 lbs.
Mattress Weight Limit: 500 lbs (36"); 600 lbs (39" and 42")
Warranty: 8 years, all components. Non-prorated
Flammability: All models conform to NFPA 101 (Life Safety Code™); CPSC 16 CFR Parts 1632 and 1633; Cal TB 117 & 129

¹ May be appropriate for multiple stage II or single stage III ulcers, based on full assessment of skin status and repositioning required, according to best clinical practice and judgment

+ Hill-Rom is a mark, and CareAssist® is a registered mark, of Hill-Rom® Services, Inc.