

# Geo-Mattress® HC

## Your Best Bet in Group I Surfaces

Specifically designed to aid in pressure ulcer prevention in the home care setting. Geo-Mattress HC is cut from a single layer of high-quality, combustion modified foam in a fluid-proof, nylon and vinyl cover. Cross-cut top surface is segmented in three distinct zones. Folded and boxed to save freight and storage.


Note: This product does not include Heel Slope® feature.

Item #	Description
SP812-29	Geo-Mattress® HC, 75"L x 35"W x 5"H
SP810-29	Geo-Mattress® HC, 80"L x 35"W x 5"H
SP811-29	Geo-Mattress® HC, 84"L x 35"W x 5"H
<b>Medicare Approved:</b> E0184 (Group I)	
<b>Weight limit:</b> 250 lbs.	<b>Mattress Weight:</b> 17 lbs.
<b>Warranty:</b> 2 years, non-prorated [applies in Home Care use only]	

All models are fire barrier version (FBV) as indicated by "-29" suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.

### Description

- ✓ Therapeutic Home Care Foam Mattress
- ✓ Non-powered Pressure Redistribution



# Geo-Mattress® 350

## High-Value Option Meets Multiple Needs


Segmented top surface, 6" height and fluid-proof cover make the new 350 a great choice for heavier home care patients. Or, an economical yet effective option for long term care residents.

Item #	Description
67535-29	Geo-Mattress 350, 75"L x 35"W x 6"H
68035-29	Geo-Mattress 350, 80"L x 35"W x 6"H
68435-29	Geo-Mattress 350, 84"L x 35"W x 6"H
<b>Medicare Approved:</b> E0184 (Group I)	
<b>Weight limit:</b> 350 lbs.	<b>Mattress Weight:</b> 21 lbs.
<b>Warranty:</b> 5 years, (2 years prorated) [applies in all care settings]	

All models are fire barrier version (FBV) as indicated by "-29" suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.

### Description

- ✓ 6" Therapeutic Foam Mattress
- ✓ Non-powered Pressure Redistribution



# Gel 80™

## Your Best Bet in Group I Surfaces

The Gel-80™, an advanced Group 1 support surface gel/foam overlay, delivers effective pressure redistribution by encapsulating multiple gel bladders within a high-density sculpted support base. Features include a wipe-clean, fluid-impervious fabric top surface, elastic corner straps and one-piece construction. Folded and boxed for easy shipping and storage.

Item #	Description
GEL-80	Gel/foam overlay for hospital-sized twin (80x35) mattress
<b>Medicare Approved:</b> E0185	
<b>Weight limit:</b> 300 lbs.	<b>Dimensions:</b> 78"L x 34"W x 3.5"H
<b>Warranty:</b> 1 year	<b>Overlay Weight:</b> 46 lbs.

Components conform to applicable California Technical Bulletin 117 standards.

### Description

- ✓ Gel/Foam Mattress Overlay
- ✓ Non-powered Pressure Redistribution

