

Advantage

Smart technology for advanced clinical and resident positioning

- Smart technology synchronizes motors offering advanced clinical and resident positioning including auto-contour, Trendelenburg, Reverse Trendelenburg and Comfort Chair.
- Rises from 7-7/8" low to 30" high.



Q-Series
Quality. Detail. Value.

Item #	Description
QD1000B	Advantage bed
Q6663	Swivel/locking 360° bed casters
Q6615	Foot end floor lock system
Q6674	Foot end control - Advantage
RM930	Three position pivot assist rail - universal
RM950	Fixed assist rail - universal
RAD306H	Standard 1/2 head end side rail - set
RAD501	3/4 length collapsible side rails
Q6412LP	Side rail mount pendant holder
Q6414LP	Gooseneck pendant holder - side rail mount
Q6416LP	Plastic pendant holder - head/foot board mount
Q6417LP	Gooseneck pendant holder - board mount
Q5302	Patient helper, trapeze bar, chain and hand grip
Q6583-2	Trapeze bracket - dual sockets
CS82	IV pole adjustable

Height Floor to Mattress Support	7-7/8" - 30"
Deck Angles - Back Elevation	0 - 70°
Deck Angles - Knee Elevation	0 - 20°
Bed Length - Overall Length	84-3/4" / 88-3/4"
Bed Length - Sleep Surface	76" / 80"
Bed Width - Sleep Surface	35"
Weight of Bed (without accessories)	160 lbs / 73 kg
Safe Working Load	500 lbs / 227 kg

For ordering and technical support on beds, accessories and furniture, call: **800-268-8671**

