


# Span+Aids® Patient Positioners

Category: Head & Body




## Head Support

Helps disperse pressure around cranium. Cutout floats ear when in the lateral position, or occiput when in supine position.

SP230-000 Standard Density 12/cs  
50676-371 Economy Density 12/cs

### Pediatric:

50676-372 Economy Density 12/cs




## Geo-Matt® 30° Body Wedge

Segmented surface provides pressure dispersion while supporting trunk at a 30° angle. Also used to support arm or leg in side lying, or at foot of bed to prevent sliding in supine or gatched position.


(7 3/8" x 11-5/8" x 20")

50960-300\* 6/cs  
50960-301 2/cs (includes fire retardant Soft Skin®)  
50960-300-96 6/cs (includes vinyl cover)



## O.R. Head Positioner

Same as Head Support (above), but slotted to provide clearance for tubes during prone (face-down) procedures.  
SP124-000 12/cs



## 45° Body Wedge

Supports the entire length of the torso.

50630-507\* Standard 6/cs (9" x 10-1/2" x 24")  
50630-307 Economy 2/cs (8" x 9-1/2" x 21")



## Geo-Matt® Pillow

Comfortable support for head and neck. Disperses pressure via hundreds of segmented cells. Comes with Soft Skin® protective cover.  
50960-200 6/cs



## Wedge


For lateral positioning and turning.  
SP039-000 (16" x 17" x 9") 6/cs



## Positioning Roll

Useful for various positioning needs, including cervical or popliteal support.

SP266-000 5" Diameter, 17" L 9/cs  
50676-377 5" Diameter, 12" L 12/cs



## Wedge

For lateral positioning and turning.  
SP265-000 (18" x 8" x 8") 8/cs

\*Indicates items vacuum packed to save storage.  
All items conform to California Technical bulletin 117.

