

PressureGuard® Air Therapy Surfaces Matrix

	Non-powered Air Therapy			Alternating/Lateral Rotation			Microclimate Management ["Low Air Loss" (LAL)]			Bariatric			
	Custom Care®	Custom Care® Convertible	Custom Care® Convertible LAL	APM	APM™	Safety™ Supreme	Custom Care® Convertible w/ control unit	Protocol®	Easy Air® EasyAir® LR	Custom Care® Convertible LAL w/ control unit	Bariatric CFT®	Bariatric APM	Easy Air® XL
Therapy Modalities	✓	✓	✓								✓		
Non-powered pressure redistribution	✓	✓	✓										
Convertible to powered via add-on control unit		✓	✓										
Frame Types				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Alternating Pressure				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Lateral Rotation (P) = Programmable				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Microclimate Management (LAL)				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Powered Flotation				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Versions to accommodate CareAssist® and recessed deck frames	✓	✓	✓				✓						
Alarm: (V)visual, (A) Audible		V, A	V	V, A	V, A	V, A	V, A	V, A	V	V	V, A	V, A	V
Timed Auto Firm (20 minutes)				✓	✓	✓	✓						
Microclimate Management	P	P	P	P	P	P	P	A	A	A	P	P	A
(P) = Passive, (A) = Active				✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Geo-Matt® top surface	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Safety Edge™ bolster design (RP) = Raised perimeter	✓	✓	✓	✓	✓	✓ (RP)	✓	✓	✓	✓	✓	✓	✓
Can maintain air in power outage; (A) = Automatically	✓ (A)	✓ (A)	✓ (A)	✓	✓	✓	✓ (A)	✓	✓ (A)	✓ (A)	✓	✓	✓ (A)
Heel Slope™; (F) = Foam-only heel section	✓ (F)	✓ (F)	✓ (F)	✓	✓	✓	✓ (F)	✓	✓	✓ (F)	✓	✓	✓
100% Latex free, 100% DHEP-free	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Non-stretch nylon with butyl/neoprene coating	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Bi-directional stretch with Shear Transfer Zone™ design	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Zippered air flow coverlet with Air Diffusion Matrix™ layer								✓	✓	✓	✓	✓	✓
Cover													
Medicare Code	E0373	E0373	E0373	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0373	E0277	E0277
Therapeutic Weight Limit (lbs.)	500	500	500	350	350	350	350	400	500 (350 LR)	500 (350 LR)	750	750	1000
Warranty: Mattress, non-pro-rated	5 yrs.	5 yrs.	5 yrs.	18 mos.	18 mos.	18 mos.	5 yrs.	18 mos.	2 yrs.	5 yrs.	18 mos.	18 mos.	2 yrs.
Warranty: Control unit, non-pro-rated				18 mos.	18 mos.	18 mos.	2 yrs.	18 mos.	2 yrs.	2 yrs.	18 mos.	18 mos.	2 yrs.
Meets Fire Safety standards Cal 129, CPSC16CFRT633, Boston IX-11	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
General													

*Hill-Rom is a mark, and CareAssist® is a registered mark, of Hill-Rom® Services, Inc.