














	LOW RISK • Braden 15-18 • Walks or relieves pressure regularly • Sits intermittently	MODERATE RISK • Braden 13-14 • Infrequently relieves pressure • Sits moderate to long periods • Questionable sensation	HIGH RISK • Braden ≤ 12 • Poor mobility to relieve pressure • Sits for a long period of time • Insensate	TREATMENT • At highest risk (Braden ≤ 12) • Existing ulcer • History of sitting ulcer • Insensate
	Prevention Without Positioning		Prevention With Positioning	
Isch-Dish® 				Patented open pocket design, fitted to the user's ischial span.
Equalizer® 			Viscoelastic top, soft polymer inserts, extra-soft ischial/ coccyx relief zone, shaped foundation.	
Isch-Dish® Thin 			Low profile version of standard Isch-Dish for users requiring reduced height.	
EZ Dish® 		Simple, non-fitted design provides partial offloading of ischials & coccyx.		
Geo-Matt® Contour 		Supportive, contoured base with premium, anti-shearing Geo-Matt top layer.		
Geo-Matt® PRT 		Low profile (3"), premium density version of standard Geo-Matt cushion.	Seat/Back Combination Products For Prevention with Positioning Low to Moderate Risk Short-WaveSM Seat and Back Cushion  Geo-WaveTM Specialty Chair Cushion 	
Geo-Matt® Wedge 		Geo-Matt surface, with 4" (front) to 2" (back) taper. helps keep user's hips from sliding forward.		
Gel-T® 		Segmented gel-filled bladder cradled between 2 soft foam layers.		
Geo-Matt® Cushion 	Unique geometry creates individually responsive, anti-shearing cells in a 4" standard density cushion.			