

Span F500P

Full Body Patient Lift



For users
Requiring
**MAXIMUM
ASSISTANCE**

- **High Capacity:** 500 lb (228 kg) safe working load.
- **Care Flexibility:** Can be used to lift patient from floor; fits beneath most low beds.
- **Safe and Reliable:** Redundant controls; overload protection; emergency lowering; 3-year frame warranty.
- **Easy operation:** Simple charging; audible/visual indicators; locking rear casters.

SPAN
better science. better care.



LIFTS



SURFACES



BEDS



SEATING



FALL PROTECTION



POSITIONERS



SKIN CARE



F500P Full Body Patient Lift

Designed to be the workhorse of your Safe Patient Handling and Mobility program. Delivers a high-value combination of patient safety, caregiver protection, ease of operation and low cost of ownership.



6-point spreader bar maximizes sling options



LCD battery level display; audible low-battery indicator



2-button pendant with LED indicator



Underbed clearance, allows use with most low beds

SPECIFICATIONS	
Weight Capacity	500 lb / 228 kg
Base Legs Length	44.9"
Base Closed External Width	24.9"
Base Open External Width	42.9"
Mast Height	48.3"
Base Height	4.5"
Min. lifting height	25.8"
Max. lifting height	69.1"
Net Weight	99 lb / 45 kg
Shipping Weight	108 lb / 49 kg

- **Easy base width adjustment** with locked position; accommodates chairs and commodes.
- **Redundant operation** via control unit, allowing use if pendant is lost or damaged.
- **Floor extraction:** Can be used to lift patients from floor.
- **Overload protection** with visual/audible alert, halts lifting when capacity is exceeded.
- **Manual emergency lowering** for added safety.
- **Quick-release battery** can be charged via charging cradle (optional), or via plug-in charger (included) for charging in place.
- **Backlit LCD display** during operation for ease of viewing.
- **Simple, fast assembly** using tools provided.
- **Certified to:** ISO 10535; IEC/EN/ANSI/AAMI/CAN/CSA-C22.2 60601-1, 60101-1-2, and 60601-1-11; SGS Mark for Safety.

Item #	Description
F500P	Full Body Patient Lift
Warranty:	3 years (frame and spreader bar); 2 years (electronics); 1 year (battery)

Part of Span's line of Safe Patient Handling and Mobility solutions, including:

F500S Sit-to-Stand Lift



Floor Lifts

F450T Standing Transfer Aid



Fixed Lithium-ion



Overhead Lifts

Portable Lithium-ion

