



In both contemporary and traditional models, Kingston motion chairs deliver the ultimate in comfort and performance.

KINGSTON MOTION CHAIRS

Kingston Motion Chairs




Ottoman • Model 1351



Kingston Contemporary
Model 1385

Kingston Motion Chairs Product Specifications

Kingston
traditional or contemporary
motion chairs in a variety
of colours and fabrics.



Model: **Traditional 1375** • **Contemporary 1385**

Chair Width	24 ³ / ₄ "	24 ³ / ₄ "
Chair Depth	32"	32"
Chair Height	45"	45"
Seat Width	21"	21"
Seat Depth	19"	19"
Seat Height	18 ¹ / ₂ "	18 ¹ / ₂ "
Arm Height	24 ³ / ₄ "	24 ³ / ₄ "

Specifications and dimensions are approximate and subject to change without notice.

M.C. HEALTHCARE PRODUCTS INC.

4658 Ontario Street, Beamsville, Ontario, Canada L0R 1B4
Telephone (905) 563-8264 • Toll Free 1-800-268-8671
Facsimile (905) 563-8680 • Email beds@mchealthcare.com