

Span+Aids® Patient Positioners

Category: Extremities



Foot Drop Stop – Standard

Elevates heels, supports ankle position. Adjustable strap can accommodate ankle contracture, protect against foot drop. Square base provides stability. Ring-of-Air® design provides ventilation.
50638-511 *6/cs



Limb Elevator

Elevates casted or uncasted leg. Can be used in both supine and side lying positions.
50632-508 *6/cs
50632-308 2/cs




Foot Drop Stop – Original

Similar to standard version (above), except hexagonal base maximizes comfort by allowing user to lie with leg in multiple positions.
55555-01 *6/cs
50638-311 6/cs




Flared Leg Extension

Elevates leg and knee in extended position.
SP291-000 2/cs



Cradle Boot – Original

Floats heel, supporting foot in neutral position, prevents external leg rotation.
55555-02 *6/cs



Leg Trough

Immobilizes leg and knee in neutral position. Base is contoured to posterior leg.
SP033-000 2/cs



Cradle Boot – Standard

Floats heel, supporting foot in the neutral position. Well-ventilated boot is easy to use and apply. Square base provides stability, discouraging leg rotation.
50644-515 *6/cs



Upper Limb Positioner

Positions the entire upper arm with sloped groove for stability and comfort. Initially designed for use with burn patients.
SP121-000 2/cs



Cock-Up Splint

Designed to help maintain hand and wrist in a desired position. Holds fingers and thumb opposed to prevent contractures. Hand stirrup is adjustable.
50640-512 *8/cs

Span+Aids® Patient Positioners

Category: Surgical & Post-Operative

NOTE: In addition to the specialty surgical positioners listed below, many other Span Positioners can also be used for site-specific positioning during surgical procedures.



Open Heart Surgery Kit

Holds each leg in externally rotated position for vein presentation. Special cuts allow for drainage tubes and sensors.

- 50666-525 *Open Heart Kit (each kit includes 2 leg positioners, 2 heel protectors, 2 leg toppers) 2 kits/cs
- 50666-625 *Replacement toppers 2 pairs/cs
- 50676-376 Economy Kit (2 leg positioners only)



O.R. Table Pillow

Used during surgery in supine or lateral position to provide support for sensitive areas of the head. Segmented squares offer maximum pressure dispersion, especially during prolonged surgical procedures.

50654-519 10/cs

(other head products page 34)



Total Knee Leg Elevator

Used in knee replacement surgery, elevates leg in flexed position.

SP195-000 1/cs



Axilla Cuff

Supports and protects axilla and upper arm during surgery.

SP010-000 12/cs



90° Leg Elevator

Bends knee at 90° angle. For use during surgery or physical therapy.

SP196-000 1/cs



Thoracotomy Arm Sling – Standard

Used with an ether screen bar, it positions hand and wrist for access to thoracic cavity without occluding IV. Adjustable hand strap allows multiple wrist positions.

50668-526 *12/cs



Laminectomy Arm Cradle Set

For use in prone position. Restricts rotation of arm, cushions axilla and provides firm support for entire limb.

50670-527 *12 pairs/cs




Thoracotomy Arm Sling – Economy

(Does not include hand strap.)

50676-373 12/cs

*Indicates items vacuum packed to save storage. All items conform to California Technical bulletin 117.



Span+Aids® Patient Positioners

Category: Surgical & Post-Operative




Heel Protector

Provides comfort and protection for low to moderate risk users. Adjustable strap allows a custom fit.
SP217-000 20 pairs/cs



Ulnar Nerve Protector

A full six inches wide, the ulnar nerve protector provides protection to ulnar nerve and elbow. Adjustable ventilated design.
50675-010 12 pairs/cs



O.R. Armboard Pads

Comfortably supports entire arm during surgery when used on O.R. table armboards.
50676-378 12 pairs/cs



Sacral Pad

Cushions sacral area during surgery. Segmented squares in center help disperse pressure.
SP019-000 12 pairs/cs



Elbow Protector

Comfortable protection with adjustable strap.
SP190-000 20 pairs/cs



Foam Positioning Rings

All-purpose pressure dispersion and comfort.

SP238-000	7"Wx1.75"H	48/cs
50676-374	8"Wx3"H	12/cs
SP239-000	9"Wx1.75"H	32/cs
SP167-000	8"Wx1.75"H	24/cs

Category: Abduction Pillows



Standard Abduction Pillow

Adjustable Velfoam straps.

50642-520	*SML	(17x12x4x7)	6/cs
50642-540	*MED	(22x14x4x7)	6/cs
50642-560	*LRG	(25x16x7x7)	6/cs

Economy Abduction Pillow

Adjustable Velfoam straps.

50642-320	SML	(17x12x4x7)	2/cs
50642-340	MED	(22x14x5x7)	2/cs
50642-360	LRG	(25x16x7x7)	2/cs



Abduction Pillow with Contour & Comfort Band

Adjustable Velfoam straps, contoured leg channels and a cotton comfort band.

SP122-000	SML	(18x12x4x7)	6/cs
SP119-000	MED	(23x14x4x7)	6/cs
SP197-000	WIDE	(20x19x7x7)	2/cs
SP131-000	LRG	(25x16x6x7)	6/cs



Abduction Pillow with Comfort Band

Adjustable Velfoam straps with a cotton comfort band.

50642-820	*SML	(17x12x4x7)	6/cs
50642-840	*MED	(22x15x5x7)	6/cs
50642-860	*LRG	(25x16x7x7)	6/cs



Abduction Pillow with Contour


Adjustable Velfoam straps and contoured leg channels.

SP357-000	*SML	(17x12x4x7)	6/cs
SP358-000	*MED	(22x15x5x7)	6/cs
SP359-000	*LRG	(25x16x7x7)	6/cs
SP204-000	SML	(17x12x4x7)	2/cs
SP205-000	MED	(22x15x5x7)	2/cs
SP206-000	LRG	(25x16x7x7)	2/cs

Size Guide


If patient is:
less than 5' tall use SML
5'1" to 5'11" use MED
6'0" to 6'5" use LRG

*Indicates items vacuum packed to save storage.
All items conform to California Technical bulletin 117.



Span+Aids® Patient Positioners

Category: Head & Body




Head Support

Helps disperse pressure around cranium. Cutout floats ear when in the lateral position, or occiput when in supine position.

SP230-000 Standard Density 12/cs
50676-371 Economy Density 12/cs

Pediatric:

50676-372 Economy Density 12/cs



Geo-Matt® 30° Body Wedge

Segmented surface provides pressure dispersion while supporting trunk at a 30° angle. Also used to support arm or leg in side lying, or at foot of bed to prevent sliding in supine or gatched position.

(7 3/8" x 11-5/8" x 20")


50960-300* 6/cs
50960-301 2/cs (includes fire retardant Soft Skin®)
50960-300-96 6/cs (includes vinyl cover)



O.R. Head Positioner

Same as Head Support (above), but slotted to provide clearance for tubes during prone (face-down) procedures.

SP124-000 12/cs



45° Body Wedge

Supports the entire length of the torso.

50630-507* Standard 6/cs (9" x 10-1/2" x 24")
50630-307 Economy 2/cs (8" x 9-1/2" x 21")



Geo-Matt® Pillow

Comfortable support for head and neck. Disperses pressure via hundreds of segmented cells. Comes with Soft Skin® protective cover.

50960-200 6/cs



Wedge

For lateral positioning and turning.

SP039-000 (16" x 17" x 9") 6/cs



Positioning Roll

Useful for various positioning needs, including cervical or popliteal support.

SP266-000 5" Diameter, 17" L 9/cs
50676-377 5" Diameter, 12" L 12/cs



Wedge

For lateral positioning and turning.

SP265-000 (18" x 8" x 8") 8/cs

*Indicates items vacuum packed to save storage.
All items conform to California Technical bulletin 117.

