



RECLINERS-LOUNGE SEATING-BENCHES



DR 8183 Chippendale Recliner



DR 0383 Queen Anne Recliner



DR 0024 Small Wing Chair



DR 0044 Lounge Chair
DR 0544 Side Chair



DR 304 Chippendale Bench



DR 404 Queen Anne Bench



DR 604 Queen Anne Bench



DR 0387 Lounge Chair



DR 7104 Martha Washington Lounge Chair
DR 1204 Martha Washington Lounge Chair
with clean out



DR 8104 Martha Washington Lounge Chair
DR 0204 Martha Washington Lounge Chair
with clean out

Model	Width	Depth	Seat Height	Height	Arm Height
DR 8185 Chippendale Recliner	30"	35"	19"	42"	24"
DR 0385 Queen Anne Recliner	30"	35"	19"	42"	24"
DR 304 Chippendale Bench	40"	18"	18"	-	-
DR 404 Queen Anne Bench	48"	18"	18"	-	21"
DR 604 Queen Anne Bench	43"	17"	18"	-	-
DR 0024 Small Wing Chair	27"	29"	18"	38 ¹ / ₂ "	25 ¹ / ₄ "
DR 0044 Gooseneck Arm Chair	25"	25"	18 ¹ / ₂ "	43"	26"
DR 0544 Side Chair	25"	25"	18 ¹ / ₂ "	43"	-
DR 0387 Lounge Chair	27"	29"	20"	36"	24 ¹ / ₂ "
DR 7104 Martha Washington Lounge Chair	23"	27"	18 ¹ / ₂ "	39"	25"
DR 1204 Martha Washington Lounge Chair with clean out	23"	27"	18 ¹ / ₂ "	39"	25"
DR 8104 Martha Washington Lounge Chair	24"	27"	19"	39"	27"
DR 0204 Martha Washington Lounge Chair with clean out	24"	27"	19"	39"	27"

All Recliners, Lounge Seating and Benches are available in a variety of colours and fabrics.

Dimensions are approximate and subject to change without notice.

M.C. HEALTHCARE PRODUCTS INC.

4658 Ontario Street, Beamsville, Ontario, Canada L0R 1B4

Telephone (905) 565-8264 • Toll Free 1-800-268-8671

Facsimile (905) 565-8680 • Email beds@mchealthcare.com