

Geo-Mattress with Wings[®]

Category: Therapeutic Raised Perimeter Mattress

We Added Safety — Without Compromise!



The supportive, raised perimeter of our popular Geo-Mattress with Wings can help prevent resident injury. Instead of the same exposed, hard bolsters others use, its raised perimeter features a tissue-friendly Geo-Matt surface from edge to edge. By gently urging the user away from the edge of the bed, its tissue-friendly raised perimeter can help reduce the risk of accidental falls. And when bed rails are in place, Geo-Mattress with Wings can help keep the user away from potentially hazardous positions between the rail and the mattress, as well as from contact with the rails themselves.

Boasting all the key benefits of the *Max*, *Plus*, and *Pro*, Geo-Mattress with Wings means safety for your residents — and their skin!



Lowered, padded openings make transfers easier and safer for user and caregiver.



Stable, high-density angled bolster, topped with patented Geo-Matt segmented surface.

QUICK FACTS

Geo-Mattress[®] with Wings[™]

Primary Use:

- Minimize fall risk due to independent exits from bed

Patient Indications:

- General patient population
- Moderate – High Risk for skin breakdown (e.g., Braden score ≥ 10)
- Decreased mobility
- Repositions frequently by self or caregiver

Key Advantages:

- Edge-to-edge antishearing surface
- Extra-firm edge bolster
- Shearing protection at transfer opening
- Heel Slope[®]

Item #	Description
W7535-29	Geo-Mattress with Wings, 75" L x 35" W x 6" H center, 8" H side
W8035-29	Geo-Mattress with Wings, 80" L x 35" W x 6" H center, 8" H side
W8435-29	Geo-Mattress with Wings, 84" L x 35" W x 6" H center, 8" H side
Weight limit:	400 lbs.
Warranty:	5 years, non-prorated
Mattress Weight:	24 lbs.

All models are fire barrier version (FBV) as indicated by "-29" suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.