

Geo-Mattress® HC

Category: Therapeutic Home Care Foam Mattress



Your Best Bet in Group I Surfaces

Specifically designed to aid in pressure ulcer prevention in the home care setting. Geo-Mattress HC is cut from a single layer of high-quality, combustion modified foam in a fluid-proof, nylon and vinyl cover. Cross-cut top surface is segmented in three distinct zones. Folded and boxed to save freight and storage.

Note: This product does not include Heel Slope® feature.


Item #	Description
SP812-29	Geo-Mattress® HC, 75"L x 35"W x 5"H
SP810-29	Geo-Mattress® HC, 80"L x 35"W x 5"H
SP811-29	Geo-Mattress® HC, 84"L x 35"W x 5"H

Medicare Approved: E0184 (Group I)
Weight limit: 250 lbs.
Warranty: 2 years, non-prorated
 [applies in Home Care use only]
Mattress Weight: 17 lbs.

All models are fire barrier version (FBV) as indicated by "-29" suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.

Geo-Mattress® 350

Category: 6" Therapeutic Foam Mattress



High-Value Option Meets Multiple Needs

Segmented top surface, 6" height and fluid-proof cover make the new 350 a great choice for heavier home care patients. Or, an economical yet effective option for long term care residents.

Item #	Description
67535-29	Geo-Mattress 350, 75"L x 35"W x 6"H
68035-29	Geo-Mattress 350, 80"L x 35"W x 6"H
68435-29	Geo-Mattress 350, 84"L x 35"W x 6"H

Medicare Approved: E0184 (Group I)
Weight limit: 350 lbs.
Warranty: 5 years, (2 years prorated)
 [applies in all care settings]
Mattress Weight: 21 lbs.

All models are fire barrier version (FBV) as indicated by "-29" suffix, denoting certified conformance to NFPA 101 (Life Safety Code); Cal Tech Bulletin #117; 16 CFR Parts 1632 and 1633; and BFD IX-11.